NU-GEL™

HYDROGEL WITH ALGINATE

HYDROGEL WITH ALGINATE

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IDROGEL CON ALGINATO

HIDROGEL CON ALGINATO

HYDROGEL MED ALGINAT

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HIDROGEL COM ALGINATO

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


ALGINATLI HİDROJEL

EN Product Description

NU-GEL™ Hydrogel is a transparent hydro active amorphous gel containing sodium alginate. The hydrogel creates a moist wound healing environment which assists with natural autolytic debridement whilst the alginate component enhances its absorptive capabilities. The gel can be used to soften and hydrate eschar by facilitating rehydration of the wound.

Indications

NU-GEL™ Hydrogel is indicated for debridement and desloughing of wounds together with the management of chronic wounds throughout all stages of the healing process. [Deviation for Spanish: NU-GEL™ Hydrogel is indicated for the treatment of decubitus ulcers and varicose in all states of healing.

Precautions

If wound infection develops, appropriate antimicrobial therapy should be initiated. The use of NU-GEL™ Hydrogel may be continued under medical supervision.

NU-GEL™ Hydrogel contains propylene glycol which might cause irritation or sensitisation in a small number of patients; if such reactions are noted, use of NU-GEL™ Hydrogel should be discontinued.

DIRECTIONS FOR USE

Application – see diagram, also on back of ampoule label

1. Carefully cleanse the wound.
2. Remove the label from the concertina pack by pulling the label as indicated.
3. Remove the cap, and swab the nozzle with a suitable

- antiseptic swab (A-B).
4. Twist to remove the tab (C).

NU-GEL™ Hydrogel can be applied to the wound with a simple one-handed technique (D).

- Hold the concertina pack in one hand and direct NU-GEL™ Hydrogel into the wound site by applying an even pressure to the base of the concertina pack with thumb or heel of hand. The nozzle tip should be held clear of the wound surface.
- Apply NU-GEL™ Hydrogel to the whole wound to a depth of at least 5mm.
- Cover with a suitable secondary dressing appropriate to the level of exudate.
- Any unused gel should be discarded.

Dressing Change and Removal

- Irrigate the wound with sterile saline to cleanse wound site

Depending upon the amount of wound exudate NU-GEL™ Hydrogel may be left undisturbed in the wound for up to 3 days. NU-GEL™ Hydrogel should be changed if leakage occurs through the secondary dressing.

Do not use if package is damaged.

Do not re-use (single use).

Do not re-sterilize.

The use by date of this product is printed on the packaging.

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©2018 Systagenix Wound Management Limited,