



Acelity™

TIELLE ESSENTIAL™ SILICONE BORDER

SILICONE ADHESIVE FOAM DRESSING

EN Product Description

TIELLE ESSENTIAL™ SILICONE BORDER, Silicone Adhesive Foam dressing is a sterile wound dressing, consisting of hydrophilic, absorbent, polyurethane foam, designed to absorb exudate.

The outer layer of the dressing is a pink, low friction, waterproof, polyurethane film which provides a bacterial barrier.

The perforated, silicone adhesive, wound contact layer enables exudate uptake, and prevents excess fluid causing maceration to healthy skin. The perforated wound contact layer prevents the dressing from sticking to the wound and reduces pain during dressing changes, by preventing the formation of granulation tissue into the pores of the foam, thus reducing trauma on dressing removal.

TIELLE ESSENTIAL™ SILICONE BORDER, Silicone Adhesive Foam dressing maintains a moist environment to assist in the healing process of moderate to heavily exuding wounds.

TIELLE ESSENTIAL™ SILICONE BORDER, Silicone Adhesive Foam dressing provides a gentle but secure adhesion to skin with minimal epidermal stripping and pain on removal.

Not made with natural rubber latex.

Indications

TIELLE ESSENTIAL™ SILICONE BORDER, Silicone Adhesive Foam dressing is indicated for use on moderate to heavily exuding chronic and acute wounds during the healing process. The dressing may be used on the following wounds: pressure ulcers, venous and arterial leg ulcers, diabetic ulcers, post-operative surgical wounds, trauma wounds, superficial and partial thickness burns, skin abrasions and cavity wounds as a secondary dressing.

TIELLE ESSENTIAL™ SILICONE BORDER, Silicone Adhesive Foam dressings are suitable for use under compression bandaging.

Contra-indications

TIELLE ESSENTIAL™ SILICONE BORDER, Silicone Adhesive Foam dressings are not indicated for use on the following:-

Do not use on individuals with a known sensitivity to polyurethane films, foams or silicone.

Do not use with oxidising solutions such as hypochlorite or hydrogen peroxide, as these can break down the absorbent polyurethane component of the dressing.

Surgical implantation.

Precautions

Do not use on dry or lightly exuding wounds.

Do not re-use in whole or in part, as it may compromise sterility and/or the performance of the dressing.

The dressing must not be ingested and must be kept away from children and animals.

For use under the guidance of a healthcare professional.

Directions for Use

Site Preparation

1. If necessary cleanse the wound area in accordance with standard protocols.
2. Dry the surrounding areas of the wound thoroughly before dressing application.

Dressing Application

1. Select a size of TIELLE ESSENTIAL™ SILICONE BORDER, Silicone Adhesive Foam dressing making sure the foam pad is slightly larger than

the wound. The foam pad size should extend well beyond the open wound [at least 2cm] in all directions.

2. Centre the dressing and apply gently over the wound site (white side to the wound).
3. Press the dressing on the surrounding skin by gentle application of pressure.

Dressing Change and Removal

1. TIELLE ESSENTIAL™ SILICONE BORDER, Silicone Adhesive Foam dressing can remain *in situ* up to 7 days, dependent on wound condition and the level of exudate. It may be necessary to change the dressing more frequently or whenever good wound care practice dictates.
2. Gently remove the dressing from the wound bed and dispose of according to local procedures and guidelines.
3. If required, cleanse the wound site with a suitable wound cleanser prior to application of a new dressing.

Caution: Federal Law (USA) restricts this device to sale by or on the order of a properly licensed healthcare practitioner.
This caution is not applicable outside the U.S.

Do not use if package is damaged.

Do not re-use.

Do not re-sterilize.

The use by date of this product is printed on the packaging.

SYMBOLS USED ON LABELLING



= Do not use if package is damaged.



= Do not re-use (single use).



= Do not re-sterilize.



= Sterilized using ethylene oxide.



= Batch number.



= Use-by: year and month /
Use-by: year, month, day.



= Consult instructions for use.



= Store below 25°C (77°F).



= Caution.



= Catalogue number.



= See instructions for use for exudate level.



= Keep dry.



= Depicting dressing size and absorbent island size at the widest point.

10013424

Distributed by:

©2016 Systagenix Wound Management Limited,
Gargrave, North Yorkshire, BD23 3RX, UK

Made in the UK

2016-09

 **Systagenix**

An Acelyty Company

