Let's Promote™

FIBRACOL® PLUS

COLLAGEN WOUND DRESSING WITH ALGINATE

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APÓSITO DE COLÁGENO CON ALGINATO PARA HERIDAS

ALJİNATLI KOLAJEN YARA SARGISI

Natural product - approximate dimensions

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 Systagenix

EN Product Description

FIBRACOL® Plus Collagen Wound Dressing with Alginate is an advanced wound care device composed of collagen and calcium alginate fibers. Its unique combination of natural biopolymers created by a patented process combines the structural support of collagen and the gel forming properties of alginates into a sterile, soft, absorbent, conformable topical wound dressing. In the presence of wound fluid FIBRACOL® Plus dressing maintains a physiologically moist microenvironment at the wound surface that is conducive to granulation tissue formation, epithelialization, and enables healing to proceed at a rapid rate. FIBRACOL® Plus dressing is versatile as a primary wound dressing, it can be cut to the exact size of the wound, multi-layered for the management of deep wounds and used in combination with either a semi-occlusive or non-occlusive secondary dressing.

Indications

FIBRACOL® Plus dressing is indicated for the management of exuding wounds including:

- Full-thickness and partial-thickness wounds
- Pressure ulcers
- Venous ulcers
- Ulcers caused by mixed vascular etiologies
- Diabetic ulcers
- Second-degree burns
- Donor sites and other bleeding surface wounds
- Abrasions
- Traumatic wounds healing by secondary intention
- Dehisced surgical incisions

Precautions

FIBRACOL® Plus dressing may be used when visible signs of infection are present in the wound area only when proper medical treatment addresses the underlying cause. FIBRACOL® Plus dressing may be used under compression therapy with healthcare professional supervision.

Contraindications

FIBRACOL® Plus dressing is not indicated for wounds with active vasculitis, third-degree burns, or patients with known sensitivity to collagen or alginates.

Adverse Reactions

FIBRACOL® Plus dressing should not be used on patients with known sensitivities to collagen or alginates. Discontinue use if signs of sensitivity appear.

Directions for Use

- Debride when necessary and irrigate the wound site with normal saline solution.
- Remove excess solution from surrounding skin.
- For heavily exuding wounds, apply to wound bed directly. For wounds with minimal exudate, apply to moistened wound bed; this will initiate the gel forming process.
- Pack deep wounds loosely. The dressing can be cut to size with sterile scissors. The amount of FIBRACOL® Plus dressing to be used depends on the size of the wound and the amount of exudate.
- FIBRACOL® Plus dressing may be covered with either a non-occlusive secondary dressing and fixed to the skin with a non-irritating tape or a semi-occlusive dressing (e.g TIELLE® Hydropolymer Adhesive Dressing or TIELLE® Plus Hydropolymer Adhesive Dressing).
- Reapply FIBRACOL® Plus dressing when the secondary dressing has reached its absorbent capacity or whenever good wound care practice dictates that the dressing should be changed. A heavily exuding wound may require daily or twice daily dressing changes. More moderately exuding wounds will require less frequent changes (every 2 to 4 days or as directed by a healthcare professional).

Following the initial application, irrigate the wound with saline solution. Reapply FIBRACOL® Plus dressing as previously instructed.

Dressing Removal

After gently removing the secondary dressing, lift any FIBRACOL® Plus dressing that has not formed a gel and discard. Using normal saline, gently irrigate the wound to remove any residual gel.

Do not re-use.

Do not resterilize.

Do not use if package is damaged.

The use by date of this product is printed on the packaging.

Caution: Federal Law (USA) restricts this device to sale by or on the order of a properly licensed healthcare practitioner.

This caution is not applicable outside the U.S.

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©2018 Systagenix Wound Management Limited,