BIOSORB™

GELLING FIBRE DRESSING

PANSEMENT FIBRE GÉLIFIANT

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GEL-VEZELVERBAND

PENSO DE FIBRA GELIFICANTE


GEELIKUITUSIDOS

GELERENDE FIBERFORBINDING

GELLING-FIBERFORBINDING

Acelity™

EN Product Description

BIOSORB™ Gelling Fibre Dressing is a soft, conformable non-woven dressing made from sodium carboxymethyl cellulose and strengthening cellulose fibre(s). When it comes into contact with wound exudate, the absorbent dressing forms a gel which helps to create a moist wound environment. This environment is conducive to autolytic debridement and supporting the healing process. Through the gel formation, any wound debris that is retained inside the fibre dressing, may be removed when the dressing is changed. Even when the dressing is moist, as a gelled fibre, its structure remains intact. The high vertical absorption of exudate into the fibre dressing protects the wound edge,¹ thus supporting the healing process.

Indications

BIOSORB™ Gelling Fibre Dressing is intended for use in the management of moderate to heavily exuding acute or chronic wounds including:

- Lower leg ulcers, pressure ulcers (Stage II to IV) and diabetic ulcers.
- Surgical wounds (e.g. post-operative, wounds left to heal by secondary intent and donor sites).
- Partial thickness burns.
- Traumatic wounds (e.g. abrasions and lacerations).
- Oncology wounds.

BIOSORB™ Gelling Fibre Dressing has a supportive effect in that it protects the wound edge and surrounding skin from maceration.¹

The product can also be used under compression.

When dry, BIOSORB™ Gelling Fibre Dressing can easily be cut with sterile scissors.

Contraindications

The product should not be used for people who are sensitive or allergic to the dressing and its components.

If signs of an allergic reaction are detected, discontinue the use of the product and seek advice from a healthcare professional on the next course of action.

BIOSORB™ Gelling Fibre Dressing is not intended to be used on wounds with severe bleeding, inside internal body cavities or closed wounds.

Precautions

- BIOSORB™ Gelling Fibre Dressing may be used when visible signs of infection are present in the wound area only when proper medical treatment addresses the underlying cause.
- Product is sterilized by irradiation and must not be resterilized.

- Product remains sterile unless the package is opened or damaged.
- The product is for single use only.

Directions for use

Site Preparation

- Clean the wound carefully. The skin surrounding the wound should be clean and dry.

Dressing Application

1. The size of the dressing should be selected and cut as necessary to overlap the wound edges by approximately 1cm. Keep the product away from the wound while cutting it to size.
2. Place the dressing on the moist wound, overlapping the wound edges by approximately 1cm. If there is a low amount of exudate present, moisten the dressing with sterile physiological saline solution (0.9%).
3. For deep wounds, use the ribbon variant. Pack wounds lightly, leaving approximately 2.5cm overlapping the wound edges to facilitate removal.
4. Apply a suitable secondary dressing to keep the product in place.

Dressing Change and Removal

The dressing should be changed when medically indicated (e.g. when the dressing has reached its absorbent capacity or when good wound care practice dictates a change is needed). The interval between changes should be no more than 7 days.

Should the dressing adhere to wounds with lower exudate levels, moisten it with sterile physiological saline solution (0.9%) before changing the dressing so that the healing process is not disturbed. Any gel residue on the wound should be removed when cleansing the wound. Deep wounds in particular should be well irrigated. Dispose in accordance with local guidance. Reapplication or continual use should not exceed 30 days unless under specific medical supervision.

¹in vitro data on file.

Do not use if package is damaged.

Do not re-use.

Do not resterilize.

The use by date of this product is printed on the packaging.

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©2016 Systagenix Wound Management Limited,