DYNA-FLEX®

Multi-Layer Compression System

Contents

1 ADAPTIC® Non-Adhering Dressing (sterile)

5 in x 9 in (12.7 cm x 22.9 cm)

A primary dressing made of knitted cellulose acetate fabric and impregnated with a specially formulated petrolatum emulsion. It is designed to help protect the wound while preventing the dressing from adhering to the wound.

1 DYNA-FLEX® Layer One (Padding/Absorption Layer) (non-sterile)

4 in x 156 in (10.2 cm x 396.2 cm)

1 DYNA-FLEX® Layer Two (Compression Bandage) (non-sterile)

4 in x 108 in (10.2 cm x 274.3 cm)

1 DYNA-FLEX® Layer Three (Cohesive Compression Bandage) (non-sterile)

4 in x 98 in (10.2 cm x 248.9 cm)

CAUTION: THIS PRODUCT CONTAINS NATURAL RUBBER LATEX WHICH MAY CAUSE ALLERGIC REACTIONS.

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7025

EN INDICATIONS

DYNA-FLEX® Multi-Layer Compression System is indicated for the management of venous leg ulcers and related conditions. The compression system may be used on patients with an ankle circumference 18 cm or larger (padded). Compression data indicates that the DYNA-FLEX® Multi-Layer Compression System provides sustained compression for up to seven days at both the ankle and below the knee, when used according to package insert directions.

CONTRA-INDICATIONS

The DYNA-FLEX® Multi-Layer Compression System should not be used on diabetic patients with small vessel disease, patients with lymphedema as a result of cancer, or on patients with an ankle brachial index (ABI) of less than 0.8. Doppler ultrasound is recommended to rule out arterial disease.

APPLICATION INSTRUCTIONS

1. Confirm the absence of arterial or diabetic small vessel disease. (Before first application.)
2. Confirm ankle circumference is greater than 18 cm. (Ankle circumference may change as edema is reduced.)
3. Wound contact dressing:
Carefully clean the wound following facility protocol and ensure skin surrounding the wound is dry.
Cover the wound with ADAPTIC® Non-Adhering Dressing.
NOTE: DYNA-FLEX® Multi-Layer Compression System may also be used over TIELLE® Hydropolymer Dressing, PROMOGRAN® Matrix Wound Dressing or PROMOGRAN PRISMA® Matrix.
4. DYNA-FLEX® Layer One (Padding/Absorption Layer, Blue Labeled Bandage):
To secure in place, wrap the padding/absorption layer against the skin. Use a spiral technique at 50% overlap from the base of the toes to just below the knee. Stretch this layer around bony prominence, such as the heel, to increase conformability. Cut off any excess bandage to avoid overlap at the knee. Secure end of Layer One with appropriate tape. (See Figures 1 and 2.)
5. DYNA-FLEX® Layer Two (Compression Bandage, Green Labeled Bandage):
To secure in place, wrap the compression bandage from the base of the toes to the ankle using the spiral technique. Begin using the figure eight technique at 50% extension from the ankle to just below the knee. Cut off any excess bandage to avoid overlap at the knee. To assist in wrapping it is recommended that the rectangular pattern on the bandage be stretched to a square pattern. Secure end of Layer Two with appropriate tape. (See Figure 3 and 4.)
6. DYNA-FLEX® Layer Three (Cohesive Compression Bandage, Red Labeled Bandage):
Wrap the cohesive compression bandage from the base of the toes to just below the knee using the spiral technique at 50% extension and 50% overlap. Cut off any excess bandage to avoid overlap at the knee. (See Figures 1 and 2.)
7. Removal: DYNA-FLEX® Compression System should be removed by a clinician with blunt scissors. After removal, clean and assess wound, chart progress, and continue treatment as necessary. Dispose of the contaminated dressing and wraps using the resealable ziplock package provided with the system.

SYMBOLS USED ON ADAPTIC® NON-ADHERING DRESSING

= Do not use if package is damaged.

= CE-mark and identification number of notified body.

= Do not re-use / single use.

= Batch number.

= Do not re-sterilize.

= Use by: year and month.

= Sterile. Method of sterilisation: irradiation.

= Store below 40°C / 104°F.

= Consult instructions for use.

= Manufacturer.

CAUTION: THIS PRODUCT CONTAINS NATURAL RUBBER LATEX WHICH MAY CAUSE ALLERGIC REACTIONS.

Caution: Federal Law (USA) restricts this device to sale by or on the order of a properly licensed healthcare practitioner.

This caution is not applicable outside the U.S.

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
Figure 3/Figura 3 Figure 4/Figura 4
(For layer 2)
(Para la 2.ª capa)

WARNINGS AND PRECAUTIONS

1. Users should be familiar with principles, techniques and complications of compression therapy.
2. Actual compression forces are dependent on product tension, wrapping technique and patient anatomy. If the patient has a very thin ankle or very prominent tibial crest, extra padding should be applied to these areas to prevent pressure necrosis. Clinical judgement must be used in product application.
3. Failure to detect significantly reduced arterial flow can result in pressure necrosis, ischemia and sequelae.
4. Patients should be counseled on possible outcomes and complications of compression therapy and advised of circumstances which should result in product removal or physician contact.
5. Should the patient develop pain or pale, cool or numb extremities distal to the dressing, the bandages should be promptly removed.
6. If during use of this product, it becomes wet, slips or gaps, becomes uncomfortable, creates irritation or uncharacteristic color or odor is observed from the foot, patients should consult their clinician immediately.
7. This product has components that contain natural rubber latex which may cause allergic reactions.

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Manufactured for/Fabricado para:
Systagenix Wound Management Limited,
Gargrave, North Yorkshire, BD23 3RX, U.K.

Packaged in Mexico containing products sourced from the U.S.
All contents made in the U.S. except ADAPTIC® Dressing made in the U.K.
and Layer #1 made in Belgium from EU sourced material.

US customers: to order product, or for product quality and technical questions, call 866-978-3296

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©2016 Systagenix Wound Management Limited,
Gargrave, North Yorkshire, BD23 3RX, U.K.