PROMOGRAN PRISMA™

MATRIX
MATRIZ

 **KCI™**
AN ACELITY COMPANY

PM0070K.a

EN Product description

PROMOGRAN PRISMA™ Matrix is comprised of a sterile, freeze-dried composite of 44% oxidized regenerated cellulose (ORC), 55% collagen and 1% silver-ORC. Silver-ORC contains 25 % w/w ionically bound silver, a well-known antimicrobial agent. In the presence of exudate, the PROMOGRAN PRISMA™ Matrix transforms into a soft, conformable, biodegradable gel, and thus allows contact with all areas of the wound. PROMOGRAN PRISMA™ Matrix, when covered with a semi-occlusive dressing, maintains a physiologically moist microenvironment at the wound surface. This environment is conducive to granulation tissue formation, epithelialization and optimal wound healing. PROMOGRAN PRISMA™ Matrix provides an effective antibacterial barrier as demonstrated by the in-vitro reduction of bacterial growth with common wound pathogens such as, *Pseudomonas aeruginosa*, *Staphylococcus aureus*, *Eschericia coli* and *Streptococcus pyogenes*. Reduction of bacterial bioburden in the dressing may result in reduced risk of infection. Literature reports of in-vitro testing indicate that collagen fibers provide a biodegradable matrix for cellular invasion and capillary growth. In laboratory testing, Collagen-ORC has been shown to absorb components of wound exudate.

PROMOGRAN PRISMA™ Matrix is a primary dressing that can be cut to fit wound with scissors and used in combination with either a semi-occlusive or non-occlusive secondary dressing. Prior to application in dry wounds saline solution should be used to hydrate PROMOGRAN PRISMA™ Matrix.

Indications for use

PROMOGRAN PRISMA™ Matrix is intended for the management of exuding wounds.

Under the supervision of a health care professional, PROMOGRAN PRISMA™ Matrix may be used for the management of:

- Diabetic ulcers
- Venous ulcers
- Pressure ulcers
- Ulcers caused by mixed vascular etiologies
- Full-thickness & partial thickness wounds
- Donor sites and other bleeding surface wounds
- Abrasions
- Traumatic wounds healing by secondary intention
- Dehisced surgical wounds

PROMOGRAN PRISMA™ Matrix may be used under compression therapy with healthcare professional supervision.

Precautions

PROMOGRAN PRISMA™ Matrix may be used when visible signs of infection are present in the wound area, only when proper medical treatment addresses the underlying cause. PROMOGRAN PRISMA™ Matrix is not intended to be a substitute for appropriate treatment of infection.

Clinicians/Healthcare Professionals should be aware that there is very limited data on prolonged and repeated use of silver containing dressings, particularly in children and neonates.

Contraindications

PROMOGRAN PRISMA™ Matrix is not indicated for third-degree burns, or patients with known sensitivity to silver, ORC or collagen.

DIRECTIONS FOR USE

Storage

PROMOGRAN PRISMA™ Matrix should be stored away from direct light. Over-exposure to light may cause some discolouration, however this does not affect the release of silver from the dressing. Store below 77 °F / 25 °C.

Prepare wound bed per your standard wound care protocol and debride when necessary.

For optimal effect, apply PROMOGRAN PRISMA™ Matrix directly to the whole wound bed.

- For a wound with low or no exudate apply PROMOGRAN PRISMA™ Matrix and hydrate with saline solution. This will initiate the transformation of the PROMOGRAN PRISMA™ Matrix into a gel.
- After hydration, through exposure to wound exudate or saline, the PROMOGRAN PRISMA™ Matrix gel will intimately come into contact with the wound surface.
- The biodegradable PROMOGRAN PRISMA™ Matrix gel is naturally absorbed into the body over time.
- In order to maintain a moist wound healing environment PROMOGRAN PRISMA™ Matrix must be covered with a semi-occlusive dressing (e.g. BIOCLUSIVE™ Plus Transparent Film Dressing, TIELLE™ Hydropolymer Dressing or TIELLE™ Plus Hydropolymer Dressing) or a non-occlusive secondary dressing and fixed to the skin with a non-irritating tape.
- After initial application, reapply PROMOGRAN PRISMA™ Matrix to the wound daily or per physician recommendation. It is not necessary to remove any residual PROMOGRAN PRISMA™ Matrix during dressing changes.

Do not use if package is damaged.

The use by date of this product is printed on the packaging.

Do not re-sterilize.

Caution: Federal Law (USA) restricts this device to sale by or on the order of a properly licensed healthcare practitioner.

This caution is not applicable outside the U.S.

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©2018 Systagenix Wound Management Limited,