 Let's Promote®

---

# PROMOGRAN PRISMA®

## WOUND BALANCING MATRIX

WUND-AUSGLEICHENDE MATRIX

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PM00725.g

## **EN** Product description

PROMOGRAN PRISMA® Matrix is a topically applied interactive wound therapy. The product is a sterile, freeze dried, composite of oxidised regenerated cellulose (ORC), collagen and silver-ORC, (a compound of Silver and ORC). In the presence of exudate the PROMOGRAN PRISMA® Matrix transforms into a soft and conformable, biodegradable gel, this allows contact with all areas of the wound. Saline or Ringer's solution should be used to hydrate PROMOGRAN PRISMA® Matrix on dry wounds.

PROMOGRAN PRISMA® Matrix modulates and re-balances the wound environment by the unique combination of:

1. Binding and inactivation of proteases (i.e. matrix metallo proteases, elastase and plasmin) which have been shown to be detrimental in excess in chronic wounds.
2. Binding and protection of naturally occurring growth factors against degradation by these excess proteases.

These naturally occurring, protected, growth factors are released back into the wound, while the detrimental proteases remain inactive upon bio-degradation of the PROMOGRAN PRISMA® Matrix.

Silver is a broad spectrum antimicrobial, which has been shown to be effective against wound pathogens.

PROMOGRAN PRISMA® Matrix has been shown to proliferate human dermal fibroblasts in vitro, as well as demonstrating free radical scavenging properties and anti inflammatory activity in laboratory tests.

### **Indications for use**

PROMOGRAN PRISMA® Matrix is indicated for the management of all wounds healing by secondary intent which are clear of necrotic tissue including:

- Diabetic ulcers
- Venous ulcers
- Pressure ulcers
- Ulcers caused by mixed vascular aetiologies
- Traumatic and surgical wounds

PROMOGRAN PRISMA® Matrix has known haemostatic properties.

PROMOGRAN PRISMA® Matrix can be used under compression therapy.

### **Contraindications**

PROMOGRAN PRISMA® Matrix is contraindicated in patients with known hypersensitivity to the components of this product, i.e. ORC, Collagen and Silver. Discontinue use if signs of sensitivity appear.

PROMOGRAN PRISMA® Matrix is not indicated for patients with extensive burns.

## **Caution**

Systemic antimicrobial therapy should be considered when wound infection is evident. PROMOGRAN PRISMA® Matrix may be used, under medical supervision, in conjunction with systemic antibiotics.

Clinicians/Healthcare Professionals should be aware that there are very limited data on prolonged and repeated use of silver containing dressings, particularly in children and neonates.

## **DIRECTIONS FOR USE**

### **Storage**

PROMOGRAN PRISMA® Matrix should be stored away from direct light. Over-exposure to light may cause some discolouration, however this does not impact product performance.

### **Site Preparation**

Before treatment, dry necrotic tissue must first be removed by surgical, enzymatic or autolytic debridement.

### **Application**

- For optimal effect, apply PROMOGRAN PRISMA® Matrix directly to the whole wound bed.
- For a wound with low or no exudate apply PROMOGRAN PRISMA® Matrix and hydrate with saline or Ringer's solution. This will initiate the gel forming process.
- After hydration the PROMOGRAN PRISMA® Matrix gel will intimately contact the wound.
- The biodegradable gel is naturally absorbed over time.
- PROMOGRAN PRISMA® Matrix must be covered with either gauze, a non-adhering or a hydropolymer dressing in order to maintain a moist wound healing environment.
- After initial treatment, re-treat the wound with PROMOGRAN PRISMA® Matrix up to every 72 hours depending upon the amount of exudate. It is not necessary to remove any residual PROMOGRAN PRISMA® Matrix.
- For heavily exuding or infected wounds it may be necessary to retreat the wound with PROMOGRAN PRISMA® Matrix every 24 hours.

**Do not use if package is damaged.**

**The use by date of this product is printed on the packaging.**

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©2018 Systagenix Wound Management Limited,