PREVENA PLUS™

INCISION MANAGEMENT SYSTEM

PATIENT GUIDE

KCI Negative Pressure Wound Therapy Units for Use at Home with PREVENA PLUS™ Dressings

PREVENA PLUS™ Dressings can be used with one of the KCI Therapy Units shown below:

PREVENA PLUS™ 125
Therapy Unit
(see page 5)

ACTIVA.C.™ Therapy Unit

V.A.C. SIMPLICITY™ Therapy Unit

V.A.C. FREEDOM™ Therapy Unit

For KCI V.A.C.® Therapy Unit Alerts and Alarms, see the user manual that comes with the therapy unit.

IMPORTANT: V.A.C.ULTA™, INFOV.A.C.™ and V.A.C.RX4™ Therapy Units **NOT** for Use at Home with PREVENA PLUS™ Dressings

V.A.C.ULTA™ Therapy Unit

INFOV.A.C.™ Therapy Unit

V.A.C.RX4™ Therapy Unit

You may have been prescribed a V.A.C.ULTA™, INFOV.A.C.™ or V.A.C.RX4™ Therapy Unit while you were in the hospital. **However, these units are not for home use.**

If you have been sent home with a V.A.C.ULTA™, INFOV.A.C.™ or V.A.C.RX4™ Therapy Unit, call your doctor right away and ask for one of the units shown above.

IMPORTANT: As with any prescription medical device, failure to consult a physician and carefully read and follow all therapy unit and dressing instructions prior to each use may lead to improper product performance and the potential for serious or fatal injury. Do not adjust therapy unit settings or perform therapy application without directions from/supervision by the treating physician.

This guide is not intended as a comprehensive resource. For additional information concerning proper use of the PREVENA PLUS™ Therapy System and important information including safety information, cautions and warnings, consult the instructions for use located in the device packaging.

PREVENA PLUS™ Incision Management System Patient Guide

You have been prescribed a PREVENA PLUS™ Dressing and a therapy unit to manage your incision. Please read this guide to learn more about these items and how to use them. This guide also includes important warnings for your safety. **Call your doctor or KCI at 1-800-275-4524 if you have any questions about the PREVENA PLUS™ Dressing or therapy unit.**

Do not re-use any part of the PREVENA PLUS™ Dressing or canisters. These items are for single use only. Re-use of these items could cause wound contamination, infection or failure to heal.

The PREVENA PLUS™ 125 Therapy Unit is a single use, disposable therapy device. This unit provides negative pressure at 125 mmHg Continuous. The PREVENA PLUS™ 125 Therapy Unit has a seven day lifespan and a rechargeable battery. Upon receipt of the therapy unit, the battery may not be fully charged.

Contraindication and Warnings

Contraindication

The dressing contains silver. If you are sensitive or allergic to silver, call your doctor right away.

Warnings

- Do not wrap carrying case strap or dressing tubing around neck.
- Small Parts - Choking Hazard
- The therapy unit is a prescription medical device. Please read and follow all the instructions in this guide so the product can perform properly while in use.

If you need:

- Magnetic Resonance Imaging (MRI), your doctor or nurse must remove the therapy unit. The dressing can remain in place.
- Hyperbaric Oxygen Therapy (HBO), your doctor or nurse must remove the therapy unit **and** dressing.

If defibrillation is required, a doctor or nurse will decide if the dressing should be removed.

Bleeding

Some patients may be at risk of heavy bleeding, with or without using the prescribed therapy system. This bleeding can be caused by the type of surgical procedure or by certain patient conditions.

If you see a large amount of blood near your incision or in the tubing or canister, leave the dressing in place. Turn off the therapy unit. Call your local emergency service (i.e. 911) right away.

Full Canister

The canister may fill up with fluid before your therapy is done. If you see a large amount of fluid in the canister, turn off the therapy unit. Call your doctor right away. The therapy unit will let you know when the canister is full.

Infection

The following symptoms may mean that your incision is infected. Call your doctor right away if you:

- have swelling at the dressing
- feel feverish
- have pus at the dressing
- feel an increase in soreness at the dressing
- feel itching at the dressing
- have a bad odor at the dressing
- have redness around the dressing
- feel warmth at the dressing

Allergic Reaction

The dressing may cause an allergic reaction if you are sensitive to acrylic adhesives or to silver. The following symptoms may mean you are having an allergic reaction. Call your doctor right away if you:

- notice redness
- develop a rash
- develop itching
- notice swelling
- develop hives

If you feel short of breath, your allergic reaction may be more serious. Immediately call your local emergency service (i.e. 911). Turn off the therapy unit and remove the dressing.

Contraindication and Warnings (continued)

- The dressing should only be applied or removed by or on the advice of your doctor.
- Do not modify the therapy unit or dressing. Do not connect the therapy unit or dressing to other devices that you might be using.
- The therapy unit is a medical device that monitors pressure and is not a toy. Keep away from children, pets and pests as they can damage the dressing and therapy unit and affect performance. Keep the therapy unit free of dust and lint as they can also cause damage and affect performance.
- To avoid damage to the therapy unit, use only the power supply and cord that comes with the unit.
- Keep the therapy unit away from heated surfaces.
- Equipment not suitable for use in the presence of flammable anesthetic mixture with air oxygen or nitrous oxide.
- When you receive the therapy unit, the battery may not be fully charged.

Instructions for Use

Do not turn off the therapy unit unless:

- your doctor instructs you
 - you experience bleeding
 - you have an allergic reaction
 - you need to disconnect the unit to shower or bathe
- See the **Warnings** section of this guide for more details.

Daily Activity

The therapy unit is small and portable. It can be worn beneath your clothing. Ask your doctor the following before wearing the device during daily activities.

- How active you can be.
- When you can resume your daily routine.

Sleeping

- Place the therapy unit in a position where the tubing will not become kinked or pinched.
- Ensure the therapy unit cannot be pulled off the table or dropped on the floor.

Showering and Bathing


- Do not use the therapy unit in the bath or shower.
- Do not place the therapy unit where it can fall or be pulled into a tub, shower or sink.
- Do not reach for a therapy unit that has fallen into water. Unplug the unit immediately if plugged in. Disconnect the unit from the dressing and call KCI.
- Before showering press the On / Off button to stop therapy. Close the tubing clamp and disconnect the tubing from the therapy unit.
- The clear drape is waterproof. You may wash or shower with the dressing in place. Avoid prolonged water contact with the dressing.
- When towel drying, be careful not to disrupt or damage the dressing.

Cleaning

If needed, clean the therapy unit and carrying case with a damp cloth and mild soap and water solution. Do not use bleach.


PREVENA PLUS™ 125 Therapy Unit

The PREVENA PLUS™ 125 Therapy Unit

Therapy Unit Carrying Case


Use the adjustable strap to wear the carrying case over your shoulder or across your chest. You can also wear the carrying case on your belt. Secure any excess tubing with the clips on the carrying case.

CAUTION: Do not wrap the carrying case strap, power cord or dressing tubing around your neck.

Turn Therapy On

Press and hold the On/Off button for three seconds.

Therapy Life Indicators

The therapy unit has a seven day lifespan.

CAUTION: This lifespan will start after the unit has run for one hour, non-stop. The lifespan will continue to count down even if the unit is turned off. The LEDs around the On/Off button show how many days are left in the lifespan.






Alerts

Alerts are indicated by LEDs and two beeps. The beeps may become louder, depending on the type of alert.

 **Alert Mute:** Press the Mute Button to silence the alert beeps. The alert beeps will come back on in two minutes if the alert condition has not been fixed.

- Refer to the Alert Table for instructions to fix the alert condition.
- LEDs cannot be turned off or muted (paused) by the user. Visual Alerts will only stop when the alert condition has been corrected.
- If the alert condition cannot be fixed, contact KCI at 1-800-275-4524.

Alert Table

Alert Type	How to fix
Blockage Alert 	Check for blockage:  Full Canister Blocked Tubing Closed Clamp
Leak Alert 	Check for leak:  at tubing connectors and at connection to canister. Press firmly around dressing edges. Use patch strips or excess drape to seal leak.
Low Battery Alert 	 Charging Fully Charged Power Connected (Battery icons: full green, partially green, empty yellow with red X)
System Error Alert  PREVENA PLUS™ 125 THERAPY UNIT	Turn therapy unit off then on.  Off On Continue therapy. Call KCI at 1-800-275-4524.

All LEDs flashing and two beeps every 15 seconds

Therapy Complete

After seven days of therapy, the PREVENA PLUS™ 125 Therapy Unit will shut off.

Disposal

The therapy unit, power supply, accessories and dressing should be removed by your doctor or nurse. Do not discard these items with household trash. This could violate local laws regarding hazardous waste.

Specifications

Transport and Storage Conditions:

Temperature Range:0°F (-18°C) to 140°F (60°C)

Relative Humidity Range:15 - 95%, non-condensing

Operating Conditions:

Temperature Range:41°F (5°C) to 104°F (40°C)

Relative Humidity Range:15 - 93%, non-condensing

Atmospheric Pressure Range for Optimum Performance:1060 hpa (-1253 ft/-381.9 m) to 700 hpa (9878 ft/3010 m)

The PREVENA PLUS™ Incision Management System is classified as a Type BF applied part under IEC 60601-1.

IP22 - Protection level against ingress of solid foreign objects and liquids.

All alerts are classified as low priority according to IEC 60601-1-8.

Conforms to: IEC 60601-1, IEC 60601-1-6, IEC 60601-1-11, IEC 60601-1-8

Electrical Equipment Interference

Cell phones or other similar products could affect the therapy unit. Move the therapy unit away from these devices if interference is suspected.

Explanation of Symbols Used

Do not use if package is damaged or open

IP22

Ingress Protection

Type BF applied part

STERILE R

Sterile using radiation

Date of Manufacture

Manufacturer

Fragile

Keep Dry

Each

Content Information

Class II Device

Temperature Limits

MR Unsafe

No Bathing or Showering

Rx only

CAUTION: Federal (US) law restricts this device to sale by or on the order of a physician.

Consult Instructions for Use

Refer to Clinician Guide

Single Use Only

EC REP

Authorized Representative in the European Community

Tripping Hazard

Use By

REF

Catalog Number

LOT

Lot Number

Do Not Resterilize

This product is designated for separate collection at an appropriate collection point. Do not dispose of as household waste.

Contains Phthalates (SENSAT.R.A.C.™ Pad Tubing)

Contact Information

For questions about this product, supplies, maintenance or more information about KCI products and services, please contact KCI or a KCI authorized representative, or:

In the US call 1-800-275-4524 or visit www.acelity.com

KCI USA, Inc. 12930 IH 10 West, San Antonio, Texas, 78249

Outside the US visit www.kci-medical.com

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Todas las marcas comerciales mencionadas en este documento son propiedad de KCI Licensing, Inc., sus filiales y/o licenciatarios. Copyright 2017 KCI Licensing, Inc.