 Let's Promote™

PROMOGRAN®

PROTEASE MODULATING MATRIX

PROTEASE MODULIERENDE MATRIX

MATRICE À EFFET ANTI PROTEASE

MATRICE MODULANTE DI PROTEASI

MATRIZ MODULANTE DE PROTEASAS

PROTEASMODULERANDE MATRIX

PROTEASE MODULERENDE MATRIX

MATRIZ MODULADORA DA PROTEASE

PROTEASEMODULERENDE MATRIX

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 Systagenix

EN Product description

PROMOGRAN® matrix is a topically applied interactive wound therapy. The product is a sterile, freeze dried composite of oxidised regenerated cellulose (ORC) and collagen. In the presence of exudate the PROMOGRAN® matrix transforms into a soft and conformable, biodegradable gel; this allows contact with all areas of the wound. Saline or Ringer's solution should be used to hydrate PROMOGRAN® matrix on dry wounds.

PROMOGRAN® matrix modulates and re-balances the wound environment by the unique combination of:

1. Binding and inactivation of proteases (i.e. matrix metallo proteases, elastase and plasmin) which have been shown to be detrimental in excess in chronic wounds.
2. Binding and protection of naturally occurring growth factors against degradation by these excess proteases.

These naturally occurring, protected, growth factors are released back into the wound, while the detrimental proteases remain inactive upon bio-degradation of the PROMOGRAN® matrix.

PROMOGRAN® matrix has demonstrated free radical scavenging properties.

Indications for use

PROMOGRAN® matrix is indicated for the management of all wounds healing by secondary intent which are clear of necrotic tissue, including:

- Diabetic ulcers
- Venous ulcers
- Pressure ulcers
- Ulcers caused by mixed vascular aetiologies
- Traumatic and surgical wounds

PROMOGRAN® matrix has demonstrated haemostatic properties. PROMOGRAN® matrix can be used under compression therapy.

Contraindications

PROMOGRAN® matrix is contraindicated in patients with known hypersensitivity to the components of this product, i.e. ORC and Collagen. Discontinue use if signs of sensitivity appear.

Warnings

If infection is suspected during treatment, an appropriate antimicrobial dressing such as ACTISORB® Activated Charcoal Dressing with Silver or systemic therapy should be used.

No safety issues have been raised in pressure ulcers, venous ulcers and diabetic ulcers to date.

DIRECTIONS FOR USE

Site Preparation

Before treatment, dry necrotic tissue must first be removed by surgical, enzymatic or autolytic debridement.

Application

- For optimal effect, apply PROMOGRAN® matrix directly to the whole wound bed.
- For a wound with low or no exudate apply PROMOGRAN® matrix and hydrate with saline or Ringer's solution. This will initiate the gel forming process.
- After hydration the PROMOGRAN® gel will intimately contact the wound.
- The biodegradable gel is naturally absorbed over time.
- PROMOGRAN® matrix must be covered with either gauze, a non adhering or a hydrolymer dressing in order to maintain a moist wound healing environment.
- After initial treatment, retreat the wound with PROMOGRAN® matrix up to every 72 hours depending upon the amount of exudate. It is not necessary to remove any residual PROMOGRAN® matrix.

Do not use if package is damaged.

The use by date of this product is printed on the packaging.

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©2018 Systagenix Wound Management Limited,