 Let's Comfort™

TIELLE® LITE

HYDROPOLYMER ADHESIVE DRESSING
WITH

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 Systagenix

EN Product Description

TIELLE® Lite Hydropolymer Adhesive Dressing provides an exudate handling system for wounds at the final stages of healing as well as traumatic wounds. The wound contact layer provides a non-adherent surface. The island dressing is absorbent and provides a moist environment that aids in the wound healing process and allows granulation to proceed under optimum conditions.

Indications

TIELLE® Lite Dressing is indicated for the management of low/non exuding and bleeding wounds.

TIELLE® Lite Dressing should be used under health care professional direction for the following indications:

- Superficial pressure ulcers
- Epithelialising lower extremity ulcers
 - Venous
 - Arterial
 - Mixed aetiology
- Diabetic ulcers
- Donor sites

Precautions

TIELLE® Lite Dressing is not indicated for use on the following:

- Third-degree burns
- Lesions with active vasculitis as this type of ulcer needs more frequent observations by a healthcare professional

TIELLE® Lite Dressings may be used when visible signs of infection are present in the wound area only when proper medical treatment addresses the underlying cause.

DIRECTIONS FOR USE

Site Preparation

1. Prepare the wound according to wound management protocol.
2. Ensure skin surrounding the wound is dry.

Dressing Application

The size of the dressing selected should allow the absorbent island to overlap the wound edge by approximately 1 cm.

1. Peel open the package and remove the dressing.
2. Partially peel pack side backing papers. Position absorbent island centrally over wound site and smooth in place.
3. One at a time, peel away side backing papers while smoothing adhesive border onto intact skin.

Dressing Change and Removal

It is essential for proper wound management to avoid trauma to the intact skin.

1. Change dressing when wound fluid is present at the edges of the foam pad. Do not allow exudate to accumulate under the backing. The dressing may be left in place up to 7 days depending upon the amount of exudate.
2. On removal lift one corner and carefully peel back. On fragile or friable skin, water or saline may be used to break the adhesive seal.

Do not use if package is damaged.

Do not re-use (single use).

Do not resterilize.

The use by date of this product is printed on the packaging

Caution: Federal Law (USA) restricts this device to sale by or on the order of a properly licensed healthcare practitioner.

This caution is not applicable outside the U.S.

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©2018 Systagenix Wound Management Limited,