 Let's Comfort®

TIELLE® PACKING

HYDROPOLYMER DRESSING
WITH

LIQUA  **LOCK**® Technology

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 Systagenix

EN Product Description

TIELLE® Packing Hydropolymer Dressing is a highly absorbent dressing made of polyurethane. The dressing provides a moist environment that aids in the wound healing process and allows granulation to proceed under optimum conditions. During use the dressing gently expands as it takes up exudate. It allows the wound to remain moist thus encouraging auto-debridement, which may initially increase lesion size. This is normal and to be expected prior to wound granulation.

Indications

TIELLE® Packing dressing is indicated for the management of chronic and deep exudating wounds, including the following:

- Pressure sores (all stages)
- Venous ulcers
- Diabetic foot

It is specially intended to be used for patients with

- very fragile skin surrounding wound
- deeper cavity wounds

TIELLE® Packing dressing can be used under compression bandaging.

Precautions

TIELLE® Packing dressing is not indicated for use on the following:

- Third-degree burns
- Lesions with active vasculitis
- Dressing should not be used when visible signs of infection are present, however, under proper medical treatment that addresses the cause, the dressing may be used

DIRECTIONS FOR USE

Site Preparation

Prepare the wound according to wound management protocol. (Note: Upon reapplying TIELLE® Packing dressing, if wound appears healthy and clean of necrotic debris, cleansing may not

be required between dressing changes.)

Dressing Application

1. Peel open the package.
2. Peel away top interleaf using sterile forceps.
3. Peel off the dressing from the remaining interleaf with sterile forceps and smooth it directly into the wound. Avoid overpacking, as dressing will expand during use.
4. Cover and secure with a secondary dressing of choice.

Dressing Change and Removal

It is essential for proper wound management to avoid trauma to the intact skin.

1. Change dressing when leakage of exudate occurs past the secondary dressing. The dressing may be left in place up to 7 days depending upon the amount of exudates and dressing condition.
2. To remove, gently ease dressing away from wound.
Pre-moisten the dressing with sterile saline solution, if needed.

Do not use if package is damaged.

Do not reuse.

Do not resterilize.

The use by date of this product is printed on the packaging.

Caution: Federal [U.S.] Law restricts this device to sale by or on the order of a physician or properly licensed professional].

This caution is not applicable outside the U.S.

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©2018 Systagenix Wound Management Limited,