 Let's Comfort™

---

# TIELLE® PLUS HEEL

HYDROPOLYMER ADHESIVE DRESSING  
WITH

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 Systagenix

## **EN** Product Description

TIELLE® Plus Heel Hydropolymer Adhesive Dressing is an anatomically shaped dressing for use on wounds located on the heel. It can be used also on the elbow or other parts of the body on which the anatomical shape of the dressing is appropriate. The dressing provides a dynamic exudate handling system with a highly absorbent internal wicking layer for regulating moisture from exuding wounds. During use the absorbent island gently expands as it takes up exudate. The island dressing maintains a moist environment. A moist wound environment supports the wound healing process by encouraging autolytic debridement thus enabling granulation to proceed under the optimum conditions. This may initially increase lesion size, which is normal and to be expected prior to wound granulation.

### **Indications**

TIELLE® Plus Heel Dressing is indicated for the management of exuding wounds located on the heel.

TIELLE® Plus Heel Dressing should be used under health care professional direction for the following indications:

- Pressure ulcers
- Lower extremity ulcers
  - Venous
  - Arterial
  - Mixed aetiology
- Diabetic ulcers

### **Precautions**

TIELLE® Plus Heel Dressing is not indicated for use on the following:

- Third-degree burns
- Lesions with active vasculitis as this type of ulcer needs more frequent observations by a healthcare professional.

TIELLE® Plus Heel Dressings may be used when visible signs of infection are present in the wound area together with proper medical treatment that addresses the underlying cause.

## **DIRECTIONS FOR USE**

### **Site Preparation**

1. Prepare the wound according to wound management protocol.
2. Ensure skin surrounding the wound is dry.

### **Dressing Application**

**The size of the dressing selected should allow the absorbent island to overlap the wound edge by approximately 1cm.**

1. Peel open the package and remove the dressing.
2. Peel back the central release paper. Position the absorbent island over wound site and smooth in place.
3. Remove one lateral release paper from the adhesive arm and apply the adhesive border onto intact skin prior to applying the second arm. Adjust the lateral arm of the dressing to the shape of the body area that you are dressing, overlapping the adhesive border onto the dressing if necessary.

### **Dressing Change and Removal**

It is essential for proper wound management to avoid trauma to the intact skin.

1. Change dressing when wound fluid is present at the edges of the foam pad. Do not allow exudate to accumulate under the backing. The dressing may be left in place up to 7 days depending upon the amount of exudate.
2. On removal lift one corner and carefully peel back. On fragile or friable skin, water or saline may be used to break the adhesive seal.

**Do not re-use (single use).**

**Do not resterilize.**

**Do not use if package is damaged.**

**The use by date of this product is printed on the packaging.**

**Caution:** Federal [U.S.] Law restricts this device to sale by or on the order of a physician [or properly licensed professional].

*This caution is not applicable outside the U.S.*

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©2018 Systagenix Wound Management Limited,