 Let's Comfort™

TIELLE® PLUS SACRUM

HYDROPOLYMER ADHESIVE DRESSING
WITH

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 **Systagenix**

EN Product Description

TIELLE® Plus Sacrum Hydropolymer Adhesive Dressing provides a dynamic exudate handling system with a highly absorbent internal wicking layer for regulating moisture from moderate to heavily exuding wounds. During use the absorbent island gently expands as it takes up exudate. The island dressing maintains a moist environment. A moist wound environment supports the wound healing process by encouraging autolytic debridement thus enabling granulation to proceed under the optimum conditions. This may initially increase lesion size, which is normal and to be expected prior to wound granulation.

Indications

TIELLE® Plus Sacrum Dressing is indicated for the management of moderate to heavily exuding sacral wounds and other difficult to dress areas.

TIELLE® Plus Sacrum Dressing should be used under health care professional direction for the following indications:

- Sacral pressure ulcers
- Difficult to dress areas.

Precautions

TIELLE® Plus Sacrum Dressing is not indicated for use on the following:

- Third-degree burns
- Lesions with active vasculitis as this type of ulcer needs more frequent observations by a healthcare professional.

TIELLE® Plus Sacrum Dressings may be used when visible signs of infection are present in the wound area only when proper medical treatment addresses the underlying cause.

DIRECTIONS FOR USE

Site Preparation


1. Prepare the wound according to wound management protocol.
2. Ensure skin surrounding the wound is dry.

Dressing Application


The dressing should allow the absorbent island to overlap the wound edge by approximately 1cm.

1. Peel open the package and remove the dressing.
2. Position dressing with wide adhesive border at the bottom (see diagram 1).

3. If customisation of the dressing is required, DO NOT cut beyond the mid-point of the bottom adhesive border (see diagram 2).
4. Remove centre backing paper and position over wound site ensuring the wound is covered by the hydropolymer central island. Smooth centre portion in place.
5. One at a time, peel away side backing papers while smoothing the adhesive border onto intact skin.

1

2

Dressing Change and Removal

It is essential for proper wound management to avoid trauma to the intact skin.

1. Change dressing when wound fluid is present at the edges of the foam pad. Do not allow exudate to accumulate under the backing. The dressing wear time is dependent upon the amount of exudate and patient condition.
2. On removal lift one corner and carefully peel back. On fragile or friable skin, water or saline may be used to break the adhesive seal.

Do not use if package is damaged.

Do not re-use (single use).

Do not re-sterilize.

The use by date of this product is printed on the packaging.

Caution: Federal [U.S.] Law restricts this device to sale by or on the order of a physician [or properly licensed professional].

This caution is not applicable outside the U.S.

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©2018 Systagenix Wound Management Limited,