TIELLE®

SILICONE BORDER

HYDROPOLYMER DRESSING WITH SILICONE

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Acelity™

EN Product Description

TIELLE® Silicone Border Hydropolymer Dressing with Silicone is designed to absorb exudate, help maintain a moist wound healing environment and minimise the risk of maceration in low to heavily exuding wounds through the combination of absorbency and moisture vapour transmission. The dressing consists of a perforated soft silicone gel adhesive wound contact layer, an absorbent pad consisting of hydropolymer foam and a non woven superabsorbent and a highly breathable film backing. The soft silicone wound contact layer of TIELLE® Silicone Border helps minimise pain and trauma at dressing change. The film backing is waterproof and provides a bacterial barrier. This may aid in the prevention of wound contamination.

Indications

TIELLE® Silicone Border dressing is indicated for the management of low to heavily exuding wounds.

TIELLE® Silicone Border dressing should be used under health care professional direction for the following indications:

- Lower extremity ulcers
 - Venous
 - Arterial
 - Mixed aetiology
- Diabetic ulcers
- Donor sites
- Pressure ulcers
- Post-surgical or traumatic wounds

TIELLE® Silicone Border dressing is also suitable for use under compression bandaging.

TIELLE® Silicone Border dressing is suitable for use on fragile or sensitive peri-wound skin.

Dressing Components

Polyurethane, polyacrylate, polyester, nylon and silicone.

Precautions

TIELLE® Silicone Border dressing is not indicated for use on the following:

- Third-degree burns.
- Lesions with active vasculitis

TIELLE® Silicone Border may be used when visible signs of infection are present in the wound area only when proper medical treatment addresses the underlying cause.

Do not use if the patient has a known sensitivity to any of the components of the dressing.

DIRECTIONS FOR USE

Site Preparation

Prepare the wound according to the appropriate wound management protocol. Ensure skin surrounding the wound is dry.

Dressing Application

1. Choose the appropriate size TIELLE® Silicone Border dressing to allow the absorbent island to overlap the wound edge by at least 1 cm.
2. Peel open the package and remove the dressing.
3. Partially peel back the release liners to expose the silicone wound contact layer of the dressing.
4. Position the absorbent island centrally over the wound site and whilst peeling away the release liners smooth the adhesive border onto intact skin.

TIELLE® Silicone Border helps maintain a moist wound healing environment which encourages autolytic debridement allowing granulation to proceed under optimum conditions. This may initially increase lesion size which is normal and expected prior to granulation.

Dressing Change and Removal

Dressing change frequency is dictated by good wound care practice and will depend on the condition of the wound and the surrounding skin. TIELLE® Silicone Border may be left in place for up to 7 days depending upon the amount of exudate.

Remove TIELLE® Silicone Border by taking hold of one corner and carefully peeling back to avoid trauma to the intact skin.

Do not use if package is damaged.

Do not re-use (single use).

Do not resterilize.

The use by date of this product is printed on the packaging.

Caution: Federal Law (USA) restricts this device to sale by or on the order of a properly licensed healthcare practitioner. This caution is not applicable outside the U.S.

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©2018 Systagenix Wound Management Limited,