Acelity™

EN

V.A.C.® Y-CONNECTOR

APPLICATION INSTRUCTIONS

ONLY FOR USE WITH KCI V.A.C.® THERAPY SYSTEMS

As with any prescription medical device, failure to consult a physician and carefully read and follow all therapy unit and dressing instructions and safety information prior to each use may lead to improper product performance and the potential for serious or fatal injury. Do not adjust therapy unit settings or perform therapy application without directions from or supervision by the clinical caregiver.

All disposable components of the V.A.C.[®] Therapy System are for single use only. Re-use of disposable components may result in wound contamination, infection and / or failure of the wound to heal.


The V.A.C.[®] Y-Connector is used to connect two T.R.A.C.[™] Pads or SensaT.R.A.C.[™] Pads to a single V.A.C.[®] Therapy Unit.

WARNINGS

WARNING: Indications, Contraindications, Warnings, Precautions and other important Safety Information are contained in the V.A.C.[®] Therapy Safety Information included with the V.A.C.[®] Therapy Unit.

- When multiple wound sites are being treated, T.R.A.C.[™] Technology senses only one wound site, the site connected to the Y-Connector arm with the post (Fig. 1).
- A blockage or leak on the non-post side will not be detected by the unit.
- Do not connect infected wounds with non-infected wounds through a Y-Connector.
- Do not connect wounds with different etiology in which cross contamination may occur.
- Avoid using a Y-Connector to connect wounds that would be optimally treated with differing pressure settings.
- It is not recommended to use a Y-Connector to connect grafts and / or flaps.
- Only one Y-Connector should be used for each V.A.C.[®] Therapy Unit.

Fig. 1

INSTALLATION INSTRUCTIONS

1. Remove V.A.C.[®] Canister from packaging and insert into the V.A.C.[®] Therapy Unit until it locks into place.

NOTE: *If the canister is not fully engaged, the V.A.C.[®] Therapy Unit will alarm.*

2. Connect V.A.C.[®] Canister Tubing to the V.A.C.[®] Y-Connector by twisting the tubing connector until it locks in place.

NOTE: *For specific dressing application procedures please refer to the application instructions provided in the dressing package.*

3. Connect each SensaT.R.A.C.[™] / T.R.A.C.[™] Pad Tubing to the V.A.C.[®] Y-Connector by twisting the SensaT.R.A.C.[™] / T.R.A.C.[™] Pad Tubing Connector until it locks in place.
4. Turn on power to the V.A.C.[®] Therapy Unit. The clinician should select the prescribed therapy setting.
5. Initiate V.A.C.[®] Therapy. Assess dressing to ensure seal integrity. The dressing should be collapsed. V.A.C.[®] Dressings should have a wrinkled appearance. There should be no hissing sounds. For ActiV.A.C.[®], InfoV.A.C.[®] and V.A.C. Ultra[™] Therapy Systems use the Seal Check[™] screen to verify that the rate of air leakage is below the alarm threshold. If there is any evidence of non-integrity, check SensaT.R.A.C.[™] / T.R.A.C.[™] Pad and V.A.C.[®] Drape seals, tubing connections, canister insertion, and ensure clamps are open.
6. Secure any excess tubing to prevent interference with patient mobility and to reduce the chance of patient lying on tubing.

NOTE: *If a leak source is identified, patch with additional V.A.C.[®] Drape to ensure seal integrity.*

NOTE: *If the wound is over a bony prominence or in areas where weight bearing may exert additional pressure or stress to the underlying tissues, a pressure redistribution (pressure relief) surface or device should be used to optimize patient offloading.*

7. Consider the V.A.C.* Y-Connector as an extension of canister tubing. Change the V.A.C.* Y-Connector at least once a week or more frequently, as needed, when the canister is changed. Dispose of the V.A.C.* Y-Connector, the canister tubing and the canister in accordance with specific institution protocol or state and local regulations.

EXPLANATION OF SYMBOLS USED


STERILE R Method of Sterilization - Radiation


 Do not use if package is damaged or open

 Single Use Only

 Do Not Resterilize

 Consult Instructions For Use

Rx Only **CAUTION:** Federal (US) law restricts this device to sale / rental by or on the order of a physician

 0473 Conforms with the Medical Device Directive (93/42/EEC) and has been subject to the conformity procedures laid down in the council directive

KCI USA, Inc.
San Antonio, TX
78249 USA
www.kci1.com

KCI Medical Products (UK) Ltd.
11 Nimrod Way
Wimborne, Dorset, BH21 7SH
United Kingdom
www.kci-medical.com

An **Acelity** Company

All trademarks designated herein are proprietary to KCI
Licensing, Inc. and its affiliates and / or licensors.
©2015 KCI Licensing, Inc. All rights reserved.

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Varemerker som nevnes her, tilhører KCI Licensing, Inc.,
©2015 KCI Licensing, Inc. Med enerett.