GENERAL PRECAUTIONS

- Contact your clinician if the SNAP™ Dressing leaks. If you cannot see your clinician within 24 hours, remove the SNAP™ Therapy Cartridge, SNAP™ Dressing and other dressing materials and cover the wound with sterile moist gauze or other dressing as instructed by your clinician.
- Immediately contact your clinician if you develop signs of infection such as increased swelling, fever, pain or redness at the wound area. Remove the SNAP™ Therapy Cartridge, SNAP™ Dressing, and other dressing materials and cover the wound with sterile moist gauze or other dressing as instructed by your clinician.
- Patients may shower or bathe with the SNAP™ Therapy System. If bathing, do not submerge the SNAP™ Therapy System under water.
- BIOLOCK™ Technology will gel exudate in the SNAP™ Therapy Cartridge.

ENGLISH

SNAP™ THERAPY CARTRIDGE

Patient / Caregiver Guide

EXPLANATION OF SYMBOLS USED

Consult Instructions
For Use

Manufacturer

Single Use Only

Use By

Do not use if package
is damaged or open

Keep Dry

Non-Sterile

Lot number

Date of
Manufacture

MR Unsafe

Catalog Number

Rx Only

CAUTION: Federal (US) law restricts this device to sale/rental by or on the order of a physician

EC REP

Authorized Representative in the
European Community

CE 0473

Conforms with the Medical Device Directive (93/42/EEC) and has been subject to the conformity procedures laid down in the council directive.

Rx Only

KCI™

An Acelity Company

CE 0473

KCI USA, Inc.
12930 IH 10 West
San Antonio, TX 78249 USA
1-800-275-4524
www.acelity.com

EC REP

KCI Medical Products (UK), Ltd.
11 Nimrod Way
Wimborne, Dorset
BH21 7SH
United Kingdom
www.kci-medical.com

Refer to the SNAP™ Therapy System Instructions for Use provided with each SNAP™ Dressing Kit for complete safety information and instructions for use.

This is not a complete set of instructions.

SAY HELLO TO SNAP™ THERAPY

We have designed the SNAP™ Therapy Cartridge to be easy to use, lightweight and completely silent. The device delivers negative pressure wound therapy, which may promote wound healing through the removal of small amounts of exudate, infectious material and tissue debris.

Visually inspect system at least every 8 hours.

THE SNAP™ THERAPY CARTRIDGE IS ACTIVE IF:

- **Green** capacity indicator is visible and stationary
- SNAP™ Dressing has a sucked down appearance
- SNAP™ Dressing is firm to the touch

WHEN SHOULD THE SNAP™ THERAPY CARTRIDGE BE CHANGED?

The SNAP™ Therapy Cartridge should be changed when the **red** pressure discharge indicator is visible at the top of the SNAP™ Therapy Cartridge and further treatment with the SNAP™ Therapy System is desired or indicated by your clinician.

INSTRUCTIONS FOR CHANGING THE SNAP™ THERAPY CARTRIDGE

If the **red** pressure discharge indicator is visible at the top of the SNAP™ Therapy Cartridge, the SNAP™ Therapy Cartridge will need to be replaced with a new SNAP™ Therapy Cartridge following the steps below if further treatment with the SNAP™ Therapy System is desired or indicated by your clinician:

1. Remove the SNAP™ Therapy Cartridge from the carrying clip.
2. Squeeze the release buttons on the tube fitting and pull the tube fitting out of the SNAP™ Therapy Cartridge.
3. Insert the tube fitting into the new SNAP™ Therapy Cartridge.
4. Push down the activation key to release and remove the key from the new SNAP™ Therapy Cartridge.
5. Place the SNAP™ Therapy Cartridge into carrying clip.

NOTE: When the SNAP™ Therapy Strap is placed around an extremity, take care to ensure that the SNAP™ Therapy Strap is not placed too tightly as this may cause discomfort or potentially decrease blood flow to the extremity.

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©2016 KCI Licensing, Inc. Reservados todos los derechos.