# TIELLE™

## NON ADHESIVE

NON ADHESIVE HYDROPOLYMER DRESSING

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Acelity™

## **EN** Product Description

TIELLE™ Non Adhesive hydropolymer dressing is designed to absorb exudate, help maintain a moist wound healing environment and minimise the risk of maceration in low to heavily exuding wounds through the combination of absorbency and moisture vapour transmission. The dressing consists of a hydropolymer foam, a non-woven superabsorbent layer and a highly breathable polyurethane film backing.

TIELLE™ Non Adhesive dressing can be cut to size using sterile scissors.

### **Indications**

TIELLE™ Non Adhesive dressing is indicated for the management of low to heavily exuding wounds.

TIELLE™ Non Adhesive dressing should be used under health care professional direction for the following indications:

- Lower extremity ulcers
  - Venous
  - Arterial
  - Mixed aetiology
- Diabetic ulcers
- Donor sites
- Pressure ulcers
- Post-surgical or traumatic wounds

TIELLE™ Non Adhesive dressing is also suitable for use under compression bandaging.

## **Precautions**

TIELLE™ Non Adhesive dressing is not indicated for use on the following:

- Third-degree burns
- Lesions with active vasculitis

TIELLE™ Non Adhesive dressing may be used when visible signs of infection are present in the wound area only when proper medical treatment addresses the underlying cause.

**Do not use if the patient has a known sensitivity to any of the components of the dressing.**

## **DIRECTIONS FOR USE**

### **Site Preparation**

1. Prepare the wound according to the appropriate wound management protocol.

### **Dressing Application**

1. Choose an appropriate size dressing to allow at least a 1cm border to overlap the wound edge.
2. Peel open the package and remove the dressing.
3. If required, cut the dressing to the desired size with sterile scissors before removing the release liner. Discard any remaining dressing after cutting. Do not keep dressing fragments for later use on the same patient or other patients.
4. Partially peel back the release liner to expose the wound contact layer of the dressing.
5. Position the dressing centrally over the wound site with the wound contact layer facing the wound. Fully remove and discard release liner.
6. TIELLE™ Non Adhesive dressing may be secured with an appropriate secondary dressing.
7. To avoid potential adherence as exudate levels reduce, consider using

TIELLE™ Non Adhesive dressing in combination with a non-adherent primary dressing.

TIELLE™ Non Adhesive dressing helps maintain a moist wound healing environment which encourages autolytic debridement allowing granulation to proceed under optimum conditions. This may initially increase lesion size which is normal and expected prior to granulation.

### **Dressing Change and Removal**

Dressing change frequency is dictated by good wound care practice and will depend on the condition of the wound and the surrounding skin. TIELLE™ Non Adhesive dressing may be left in place for up to 7 days depending on the amount of exudate.

Remove any secondary retention dressings. Remove TIELLE™ Non Adhesive dressing by taking hold of one corner and carefully peeling back. It is essential for proper wound management to avoid trauma to the intact skin during removal.

**Do not use if package is damaged.**

**Do not re-use.**

**Do not resterilize.**

**The use by date of this product is printed on the packaging.**

### **US only**

**Caution:** Federal Law (USA) restricts this device to sale by or on the order of a properly licensed healthcare practitioner.

This caution is not applicable outside the U.S.

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©2016 Systagenix Wound Management Limited,