



V.A.C.® THERAPY  
ESSENTIALS



# V.A.C.® Negative Pressure Wound Therapy Essentials

A Reference Guide for V.A.C.® Disposables

## V.A.C.® System Dressings

Part Numbers/  
Dressings per Case

V.A.C.® Simplace™ EX Small Dressing  
2 Spiral GranuFoam™ Dressings (7.7 x 11.2 x 1.75 cm), 2 strips of  
V.A.C.® Drape, 1 SensaT.R.A.C.™ Pad with connector, 1 disposable ruler

M8275046/5



V.A.C.® Simplace™ EX Dressing System

V.A.C.® Simplace™ EX Medium Dressing  
2 Spiral GranuFoam™ Dressings (14.7 x 17.4 x 1.75 cm), 2 drape strips and 1 sheet  
of V.A.C.® Drape, 1 SensaT.R.A.C.™ Pad with connector, 1 disposable ruler

M8275045/5



V.A.C.® Simplace™ Dressing System

V.A.C.® Simplace™ Small Dressing w/ 3M™ Tegaderm™ Drape  
2 Spiral GranuFoam™ Dressings (7.7 x 11.2 x 1.75 cm), 2 sheets of 3M™  
Tegaderm™ Drape, 1 SensaT.R.A.C.™ Pad with connector, 1 disposable ruler

M8275041/10

M8275041/5



V.A.C.® GranuFoam™  
Bridge Dressing System

V.A.C.® Simplace™ Medium Dressing w/ 3M™ Tegaderm™ Drape  
2 Spiral GranuFoam™ Dressings (14.7 x 17.4 x 1.75 cm), 3 sheets of 3M™  
Tegaderm™ Drape, 1 SensaT.R.A.C.™ Pad with connector, 1 disposable ruler

M8275040/10

M8275040/5

V.A.C.® GranuFoam™ Bridge Dressing  
1 GranuFoam™ Dressing (6 x 17 x 1.9 cm), 1 GranuFoam™ Bridge Dressing (67 cm) with  
integrated SensaT.R.A.C.™ Pad, 1 sheet of perforated V.A.C.® Drape with pre-cut hole,  
1 disposable ruler

M8275042/10

M8275042/5

V.A.C.® GranuFoam™ Bridge XG Dressing  
2 Spiral GranuFoam™ Dressings, 1 GranuFoam™ Bridge Dressing (67 cm) with integrated  
SensaT.R.A.C.™ Pad, 1 sheet of V.A.C.® Drape, 1 sheet of perforated V.A.C.® Drape with pre-  
cut hole, 1 disposable ruler

M8275044/5



V.A.C.® GranuFoam™  
Bridge XG Dressing System

V.A.C. GranuFoam Silver™ Small Dressing  
1 GranuFoam Silver™ Dressing (10 x 7.5 x 3.2 cm), 1 V.A.C.® Drape,  
1 SensaT.R.A.C.™ Pad with connector, 1 disposable ruler

M8275098/10

M8275098/5

V.A.C. GranuFoam Silver™ Medium Dressing  
1 GranuFoam Silver™ Dressing (18 x 12.5 x 3.2 cm), 2 sheets of V.A.C.® Drape,  
1 SensaT.R.A.C.™ Pad with connector, 1 disposable ruler

M8275096/10

M8275096/5



V.A.C. GranuFoam Silver™  
Dressing System

V.A.C. GranuFoam Silver™ Large Dressing  
1 GranuFoam Silver™ Dressing (26 x 15 x 3.2 cm), 2 sheets of V.A.C.® Drape,  
1 SensaT.R.A.C.™ Pad with connector, 1 disposable ruler

M8275099/10

M8275099/5

V.A.C.® GranuFoam™ Small Dressing  
1 GranuFoam™ Dressing (10 x 7.5 x 3.2 cm), 1 sheet of V.A.C.® Drape,  
1 SensaT.R.A.C.™ Pad with connector, 1 disposable ruler

M8275051/10

M8275051/5

V.A.C.® GranuFoam™ Medium Dressing  
1 GranuFoam™ Dressing (18 x 12.5 x 3.2 cm), 2 sheets of V.A.C.® Drape,  
1 SensaT.R.A.C.™ Pad with connector, 1 disposable ruler

M8275052/10

M8275052/5



V.A.C.® GranuFoam™  
Dressing System

V.A.C.® GranuFoam™ Large Dressing  
1 GranuFoam™ Dressing (26 x 15 x 3.2 cm), 2 sheets of V.A.C.® Drape,  
1 SensaT.R.A.C.™ Pad with connector, 1 disposable ruler

M8275053/10

M8275053/5

V.A.C.® GranuFoam™ X-Large Dressing  
1 GranuFoam™ Dressing (60 x 30 x 1.5 cm), 5 sheets of V.A.C.® Drape,  
1 SensaT.R.A.C.™ Pad with connector, 1 disposable ruler

M8275065/5



V.A.C.® GranuFoam™  
Round Dressing System

V.A.C.® GranuFoam™ Round Dressing  
1 perforated GranuFoam® Dressing (5" diameter), 1 sheet of V.A.C.® Drape,  
1 SensaT.R.A.C.™ Pad with connector, 1 disposable ruler

M8275075/10

M8275075/5



V.A.C.® GranuFoam™  
Thin Dressing System

V.A.C.® GranuFoam™ Thin Dressing  
1 perforated GranuFoam® Dressing (26 x 15 x 1.6 cm), 2 sheets of V.A.C.® Drape,  
1 SensaT.R.A.C.™ Pad with connector, 1 disposable ruler

M8275081/10

M8275081/5

V.A.C.® GranuFoam™ Hand Dressing  
1 hand-shaped GranuFoam® Dressing with finger separators, 1 sealing bag,  
1 SensaT.R.A.C.™ Pad with connector, 1 disposable ruler

M8275064/5



V.A.C.® GranuFoam™  
Hand Dressing System

V.A.C.® WhiteFoam Small Dressing (foam only) 1 polyvinyl alcohol WhiteFoam Dressing (10 x 7.5 x 1 cm)	M6275033/10
V.A.C.® WhiteFoam Large Dressing (foam only) 1 polyvinyl alcohol WhiteFoam Dressing (10 x 15 x 1 cm)	M6275034/10
V.A.C.® WhiteFoam Small Dressing 1 polyvinyl alcohol WhiteFoam Dressing (10 x 7.5 x 1 cm), 1 V.A.C.® Drape, 1 SensaT.R.A.C.™ Pad with connector, 1 disposable ruler	M8275068/10 M8275068/5
V.A.C.® WhiteFoam Large Dressing 1 polyvinyl alcohol WhiteFoam Dressing (10 x 15 x 1 cm), 1 V.A.C.® Drape, 1 SensaT.R.A.C.™ Pad with connector, 1 disposable ruler	M8275067/10 M8275067/5



V.A.C.® WhiteFoam Dressing



V.A.C.® WhiteFoam Dressing System

V.A.C.® System Canisters Package Contents	Part Numbers/Quantity per Case
--	--------------------------------

InfoV.A.C.® 500 mL Canister with Gel 1 canister (sterile) with SensaT.R.A.C.™ Tubing, clamp and connector (Also used with V.A.C.Ultra™ Therapy Unit)	M8275063/10 M8275063/5
InfoV.A.C.® 500 mL Canister without Gel 1 canister (sterile) with SensaT.R.A.C.™ Tubing, clamp and connector (Also used with V.A.C.Ultra™ Therapy Unit)	M8275071/10 M8275071/5
InfoV.A.C.® 1000 mL Canister with Gel* 1 canister (sterile) with SensaT.R.A.C.™ Tubing, clamp and connector (Also used with V.A.C.Ultra™ Therapy Unit)	M8275093/5
V.A.C. ATS® and V.A.C. Instill® 500 mL Canister with Gel 1 canister (sterile), tubing, clamp and connector	M6275063/10 M6275063/5
V.A.C. ATS® and V.A.C. Instill® 500 mL Canister without Gel 1 canister (sterile), tubing, clamp and connector	M6275071/10 M6275071/5
V.A.C. ATS® and V.A.C. Instill® 1000 mL Canister with Gel* 1 canister (sterile), tubing, clamp and connector	M6275093/5
ActiV.A.C.® 300 mL Canister with Gel 1 canister (sterile fluid path) with SensaT.R.A.C.™ Tubing, clamp and connector	M8275058/10 M8275058/5
V.A.C. Freedom® 300 mL Canister with Gel 1 canister (sterile fluid path), tubing, clamp and connector	320058/10 320058/5



InfoV.A.C.® 500 mL Canister



InfoV.A.C.® 1000 mL Canister



V.A.C. ATS® 500 mL Canister



V.A.C. ATS® 1000 mL Canister



ActiV.A.C.® 300 mL Canister



V.A.C. Freedom® 300 mL Canister

\*1000 mL canister is recommended for Acute Care (hospital) use only

V.A.C.® System Accessories Package Contents	Part Numbers/Quantity per Case
--	--------------------------------

SensaT.R.A.C.™ Pad Only 1 SensaT.R.A.C.™ Pad with Tubing, clamp and connector	M8275057/10
V.A.C.® Drape 1 adhesive drape (30.5 x 26 cm)	M6275009/10
V.A.C.® Y-Connector Allows two dressings to be connected to one V.A.C.® Therapy Unit	M6275066/10 M6275066/5
V.A.C.® Tubing Cap Used to create seal at end of canister tubing	M6275069/10 M6275069/5



SensaT.R.A.C.™ Pad



V.A.C.® Drape



V.A.C.® Y-Connector



V.A.C.® Tubing Cap

**ABThera™ Open Abdomen Therapy System**  
Description & Package Contents

Part Numbers/Quantity per Case

**ABThera™ SensaT.R.A.C.™ Open Abdomen Dressing**

1 Visceral Protective Layer, Perforated GranuFoam™ Dressing,  
V.A.C.® Drape, SensaT.R.A.C.™ Tubing Set

M8275026/5

ABThera™ SensaT.R.A.C.™ OA Dressing is for use with the negative pressure therapy provided by the InfoV.A.C.™ Therapy Unit and the V.A.C.Ultra™ Therapy Unit. It is not for use with V.A.C. VeraFlo™ Therapy provided by the V.A.C.Ultra™ Therapy Unit.



ABThera™ SensaT.R.A.C.™  
Open Abdomen Dressing

**Prevena™ Incision Management System**  
Description & Package Contents

Order Number

**Prevena™ Peel & Place™ Dressing, 5-Pack**

1 Prevena™ Foam Dressing with pressure indicator (10 x 3 x 1.75 cm);  
Prevena™ Patch Strips™; 1 V.A.C.® Connector

PRE1055US

Prevena™ Peel & Place™ Dressings are for use with the negative pressure provided by the InfoV.A.C.™ Therapy Unit; they are not for use with the V.A.C.Ultra™ Therapy Unit.



Prevena™ Incision Management System

**Prevena™ Customizable™ Dressing, 5-Pack**

2 Hydrocolloid Sealing Strips for each dressing; 4 V.A.C.® Drapes for each  
dressing; 1 Interface Pad for each dressing; 1 Ruler for each dressing; 1  
V.A.C.® Therapy Unit Connector for each dressing

PRE2055US

**V.A.C.Ultra™ Therapy System**  
Description & Package Contents

Order Number

The following components are for use with the V.A.C. VeraFlo™ Therapy provided  
by the V.A.C.Ultra™ Therapy Unit.

**V.A.C. VeraFlo™ Dressing, 5-Pack, Small**

2 foam dressings (11 x 8 x 1.8 cm), 2 sheets of V.A.C.® Advanced Drape, 1 foam  
quantity label, 2 packets of 3M™ Cavilon™ Barrier, 1 V.A.C. VeraT.R.A.C.™ Pad and  
Tubing Set, 1 disposable ruler

ULTVFL05SM



V.A.C. VeraFlo™ Dressing System

**V.A.C. VeraFlo™ Dressing, 5-Pack, Medium**

2 foam dressings (17 x 15 x 1.8 cm), 2 sheets of V.A.C.® Advanced Drape, 1 foam  
quantity label, 2 packets of 3M™ Cavilon™ Barrier, 1 V.A.C. VeraT.R.A.C.™ Pad and  
Tubing Set, 1 disposable ruler

ULTVFL05MD

**V.A.C. VeraFlo™ Dressing, 5-Pack, Large**

2 foam dressings (25 x 15 x 1.6 cm), 5 sheets of V.A.C.® Advanced Drape, 1  
foam quantity label, 4 packets of 3M™ Cavilon™ Barrier, 1 V.A.C. VeraT.R.A.C.  
Duo™ Pad and Tubing Set, 1 disposable ruler

ULTVFL05LG



V.A.C. VeraFlo™ Cleanse™ Dressing System

**V.A.C. VeraFlo Cleanse™ Dressing, 5-pack, Medium**

1 foam dressing (61 x 3.2 cm diameter), 2 sheets of V.A.C.® Advanced Drape,  
1 foam quantity label, 2 packets of 3M™ Cavilon™ Barrier, 1 V.A.C.  
VeraT.R.A.C.™ Pad and Tubing Set, 1 disposable ruler

ULTVCL05MD



V.A.C. VeraLink™  
Cassette

**V.A.C. VeraLink™ Cassette, 5-Pack**

1 V.A.C. VeraLink™ Cassette with instillation tubing

ULTLNK0500

**V.A.C. VeraT.R.A.C. Duo™ Tubing Set, 5-Pack**

1 V.A.C. VeraT.R.A.C. Duo™ Tubing Set

ULTDUO0500



V.A.C. VeraT.R.A.C.  
Duo™ Tubing Set

To learn more about the V.A.C.® Therapy System, please visit our website at [www.kci1.com](http://www.kci1.com) or call us at 1-800-275-4524.

Follow local institutional protocols for infection control and waste disposal procedures. Local protocols should be based on the applicable local government environmental regulations.



NOTE: Specific indications, contraindications, warnings, precautions and safety information exist for KCI products and therapies. Please consult a physician and product instructions for use prior to application. Rx only.

©2013 KCI Licensing, Inc. All rights reserved. 3M and Tegaderm are trademarks of 3M Company. All other trademarks designated herein are proprietary to KCI Licensing, Inc., its affiliates and/or licensors. DSL#15-0100.US • Lit# 29-B-113 • (REV. 8/15)

