

	Item	Product Description	Specifications	
Therapy Units	PREVENA™ 125 Therapy Unit	FDA Class II Medical Device <ul style="list-style-type: none"> Negative Pressure Options: Pre-set, continuous NPWT at -125mmHg Disposable, single patient use 3 AA batteries pre-installed 	Weight: 0.4lbs (0.20kg) Dimensions: Approx. 13.6 x 7.5 x 3cm	
	PREVENA PLUS™ 125 Therapy Unit	FDA Class II Medical Device <ul style="list-style-type: none"> Negative Pressure Options: Pre-set, continuous NPWT at -125mmHg Disposable, single patient use Rechargeable battery 	Weight with empty canister: 0.7lbs (0.32kg) Dimensions: Approx 16 x 9 x 3cm	
Dressings	<ul style="list-style-type: none"> Contain a skin friendly wicking interface layer that includes 0.019% ionic silver and polyurethane foam bolster, specially designed for application to the surgical area Silver (0.019%) in the interface layer is not intended to treat infection, but to reduce bacterial colonization in the fabric 	PEEL & PLACE™ Dressing - 20cm	Polyurethane foam bolster	25.4 x 6.4 x 1.8cm
			Polyurethane drape	35.6 x 20.3cm
			Connecting tube (integrated with dressing)	Dressing tube length: 112cm
			Appropriate incision length	Maximum 20cm (8 inches)
		PEEL & PLACE™ Dressing - 13cm	Polyurethane foam bolster	15.2 x 6.4 x 1.8cm
			Polyurethane drape	22.8 x 13.3cm
			Connecting tube (integrated with dressing)	Dressing tube length: 112cm
			Appropriate incision length	Maximum 13cm (5 inches)
		PEEL & PLACE™ Dressing - 35cm	Polyurethane foam bolster	37.8 x 6.4 x 1.8cm
			Polyurethane drape	47.5 x 13.5cm
			Dressing tube (integrated with dressing)	Dressing tube length: 112cm
			Appropriate incision length	Maximum 35cm (17 inches)
		PREVENA™ CUSTOMIZABLE™ Dressing with SENSAT.R.A.C.™ Technology	Polyurethane foam bolster with integrated hydrocolloid strips	90 x 9.8 x 1.8cm
			Hydrocolloid strips (2)	2 x 24cm each
			V.A.C.® Drape (4)	30.5 x 26cm each
			SENSAT.R.A.C.™ Pad	Tube length: 91cm
Appropriate incision length	Recommended 90cm (34.5 inches) or less			
Canisters	PREVENA™ 45ml Canister	<ul style="list-style-type: none"> 45ml canister with connecting tubing Single use, sterile, disposable 	Canister	7.2 x 6.9 x 3cm
			Canister tube	Canister tube length: 20.2cm
	PREVENA PLUS™ 150ml Canister	<ul style="list-style-type: none"> 150ml canister with connecting tubing Single use, sterile, disposable 	Canister	15.5 x 9 x 2.5cm
			PREVENA PLUS™ Connector	Tube length: 45cm
Connectors	PREVENA™ Therapy V.A.C.® Connector	Required for connection of PREVENA™ Incision Dressings to approved V.A.C.® Therapy Units	Dimensions:	2.7 x 1.5 x 1.2cm
			Weight:	1.6 grams
	PREVENA™ Y-Connector	Connects two PREVENA™ Incision Dressings to a single approved V.A.C.® Therapy Unit	Dimensions:	3.8 x 1.5 x 2.6cm
			Weight:	2.6 grams

Environmental Conditions

PREVENA™ PEEL & PLACE™ Incision Management System

Storage Temperature Range:	-4°F (-20°C) to 140°F (60°C)
Operating Temperature Range:	41°F (5°C) to 104°F (40°C)
Altitude Range for Optimum Performance:	-1253 to 9878 feet (-381.9m to 3010m)
Relative Humidity Range:	15-95%, non-condensing

PREVENA PLUS™ Incision Management System

Storage Temperature Range:	0°F (-18°C) to 140°F (60°C)
Operating Temperature Range:	41°F (5°C) to 104°F (40°C)
Altitude Range for Optimum Performance:	-1,253 to 9,878 feet (-381.9m to 3010m)
Relative Humidity Range:	15-93%, non-condensing

Electrical Specifications (PREVENA™ 125 Therapy Unit)

Battery Powered using primary (non-rechargeable) 3 AA Lithium-ion batteries

IEC Classification:

- Medical Equipment
- Electrical Safety: Internally powered
- EMC: Class B for hospital use
- Type BF Applied Part according to IEC 60601-1:2005 (3rd Ed)
- IP24: Protection against solid objects greater than 12.5mm and against liquid water falling for short periods of time
- Suitable for home use

Contains Latex	Contains DEHP (Di(2-ethylhexyl)phthalate)	Contains Mercury	Contains Haz Mat	Contains Thimerosal	Radiopaque
NO	NO	NO	NO	NO	NO

Electrical Specifications (PREVENA PLUS™ 125 Therapy Unit)

A rechargeable Lithium-ion battery pack with AC Mains disconnection IEC Classification: Medical Equipment

- Electrical Safety: Class II
- EMC: The PREVENA PLUS™ 125 Therapy Unit conforms to the intent of the 2004/108/EC EMC Directive; however, all medical equipment may produce electromagnetic interference or be susceptible to electromagnetic interference.
- Type BF Applied Part according to IEC 60601-1:2005 (3rd Ed)
- IP22: Protection against solid objects greater than 12.5mm and against liquid water falling for short periods of time
- Suitable for home use

Contains Latex	Contains DEHP (Di(2-ethylhexyl)phthalate)	Contains Mercury	Contains Haz Mat	Contains Thimerosal	Radiopaque
NO	YES, SENSAT.R.A.C.™ Pad	NO	NO	NO	NO

For more information, visit acelity.com or call **800-275-4524**

Follow local institutional protocols for infection control and waste disposal procedures. Local protocols should be based on the applicable local government environmental regulations.

NOTE: Specific indications, contraindications, warnings, precautions and safety information exist for the PREVENA™ Incision Management System. Please consult the PREVENA™ Incision Management System Clinician Guide Instructions for Use prior to application. Rx only.

Copyright 2016 KCI Licensing, Inc. All rights reserved. All trademarks designated herein are proprietary to KCI Licensing, Inc., its affiliates and/or licensors. LIT#29-A-333 • DSL#16-1242.US (08/17)